SleeperOne[®]

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Хорошо промыть и просушить
Verify that the container is free of dirt. Dry the container with a non-woven compress (or with compressed air).

Put the container in a sterile packaging as per EN 868-5.
Type B steam autoclave sterilization compliant with En13060.
Sterilisation time of minimum 18min at 275°F (134°C).

Check that the container retrieved is dry and sterile.
Store the sterile items in a dry and dustless place.

*For these different steps, follow the standard directives and specific guides of your country.

Spare parts

You can order extra accessories from your local dealer or by contacting Dental Hi TEC directly.

Plastic containers

(Ref. SA050 or SA051)

2.5m Low voltage power supply

(Ref. C2904)

Control box fixing adhesives

(Ref. SA0226)

Needles ^{DHT}

ø 0,30 (30G) - 16 mm
Ref. C3005

ø 0,30 (30G) - 9 mm
Ref. C3002

ø 0,40 (27G) - 13 mm
Ref. C3001

Technical characteristics

Name: SleeperOne Model: S4
 Device classification: IIA
 Nominal voltage: CA 100-240 V 50/60 Hz
 Ambient temperature: T_{ma}: 10-40°C
 Altitude max: 3000m
 Class II power supply enclosed: 18V - 560 mA
 Characteristics of the radio transmitter:
 Frequency / Transmission power: 868.3MHz / max. 10mW EIRP
 Data rate / Channel bandwidth / Modulation type : 120 kbps / 280 kHz / ASK

Do not use this device with a power supply other than the one supplied (FW7660).

Dimensions :
Handpiece: Diameter: 26 mm - Length: 228 mm - Weight: 70 g
Handpiece holder: Height: 75 mm - Width: 44 mm - Length: 44 mm
Control box : Length : 109 mm - Width : 78 mm - Height : 24 mm - Weight : 106 g
Wireless pedal: Length: 250 mm - Width: 180 mm - Height: 63 mm - Weight: 603g
Wire pedal (option) : Length: 185 mm - Width: 145 mm - Height: 33 mm - Weight: 320g
 Cartridge to be used: 1.8 ml cartridge (ISO 11499)
 Auto-off: After 10 minutes without use
 Compliance: 

We declare on our own liability that this product complies with the EC and CSA regulatory provisions.
 In the interest of product improvement, Dental Hi Tec reserves the right to make modifications to its products without prior notice.

SleeperOne is intended for use in an electromagnetic environment where the electromagnetic interferences are controlled and defined in tables 1 & 2 on the following pages
 The customer or user of the device must make sure that it is used in an environment as described and can help to prevent interferences by respecting minimal distance between RF communication devices (portable or mobile) (transmitters) and SleeperOne as recommended in table 3, depending on the transmission power of the communication device.

TABLE 1 - Electromagnetic emissions

Emission tests	Compliance:	Electromagnetic environment – Directives
RF emissions CISPR 11	Group 1	SleeperOne uses RF energy only for its internal function. Furthermore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
IEC 61000-3-2 Harmonics transmission	Class A	The use of SleeperOne is suitable for use in all types of establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
IEC 61000-3-3 Voltage fluctuations / Flicker emissions	Complies	

TABLE 2 - Electromagnetic immunity

Immunity test	IEC 60601- test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV ±1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV ±2 kV	
Voltage dips, short interruptions and voltage variations on power supply input lines	< 5% UT (> 95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 sec	0 V / 10 ms 92 V / 100 ms 161 V / 500 ms 0 V / 5000 ms	Mains power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during power mains interruptions, it is recommended to use an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not Applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the ac mains voltage prior to application of the test level.

Immunity test	IEC 60601 Test level required	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable or mobile RF communications equipment should not be used closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = [3,5/V1] √P d = [3,5/E1] √P 80 MHz to 800 MHz d = [7/E1] √P 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W), according to the manufacturer and d is the recommended separation distance in metres (m) ³ .
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
 a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the device should be observed to verify the normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
 b) Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

TABLE 3 - Recommended separation distances between portable & mobile RF communications equipment and SleeperOne

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 kHz to 800 MHz	800 kHz to 2,5 GHz
	$d = \left[\frac{3,5}{V1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E1} \right] \sqrt{P}$	$d = \left[\frac{7}{E1} \right] \sqrt{P}$
0,01	0,12	0,12	0,12
0,1	0,37	0,37	0,37
1	1,17	1,17	1,17
10	3,69	3,69	3,69
100	11,67	11,67	11,67

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
 NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

This appliance must not be disposed of with household waste

It contains dangerous substances, which can have important consequences for the environment and human health. It must be recycled in a specific treatment facility.

If you are concerned by the European WEEE directive, please contact your distributor or the Dental Hi Tec company who will take all of the necessary measures towards the collection and recycling of the device.

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Warranty conditions

DENTAL Hi Tec offers a two-year warranty for SleeperOne from the date of purchase.

This warranty does not apply to the sterilizable accessories. The implementation of the warranty requires the submission of the original proof of purchase indicating the date of purchase. This warranty only applies if the SleeperOne is used in compliance with the operating instructions and tolerances as indicated in this manual.

The warranty is invalid if:

- a technical repair is made to the device by an unapproved person.
- the device is damaged by an external source (lightning, fire, water damage, negligence or accident), by improper use or non-compliant installation.
- the device is adapted to operating norms that are different from those indicated on the material safety data sheet located under the control unit.
- the purchase documents are changed or altered.
- the serial number of the device has been made unreadable.

In no case shall the application of this warranty constitute an entitlement to complete replacement of the device or give the right to claim damages. SleeperOne cannot be subject to return outside of the warranty clauses.

SleeperOne®

For any complementary information,
contact our advisors:

Clinical & technical information: **+33 241 560 590**

or **sleeperone@dentalhitec.com**

DENTAL HI TEC
Innovation for dental anaesthesia

ZI de l'Appentière • 49280 Mazières-en-Mauges • FRANCE
world@dentalhitec.com

www.dentalhitec.com